



1" Turbo

A-Cure Plus

Curing Light

The new generation multiwavelength LED for producing high-intensity light

Model: A-Cure Plus

Power Input: 100-240VAC, 50Hz/60Hz

Power Output: 5VDC/2A

Wave length: 385nm-515nm

Light intensity: 800m/cm²- 2500mW/cm²

Battery voltage and capacity: 3.7V/1900mAh



A-Cure

Broad spectrum led light

Constant-current drive

Wireless & Hygienic design

Specification:

Power input 100-240VAC 50Hz/60Hz

Power output DC 5V/2A

Wave length 385nm-515nm

Light intensity 800mW/cm²-1800mW/cm²

Battery 3.7V/1900mAh

