

# VAT-3

## Activator



### Basic technical parameters

- Dimension: 165mm×28mm
- Power adapter input: AC 100-240V, 50/60Hz 0.4A
- Power adapter output: DC 5V/1A
- Main unit input: DC 3.7V 1500mAh
- Charge time: About 4 hours
- Power gear: 2 gears
- Tip vibration frequency: 45kHz±5kHz
- Tip main amplitude offset: <200um

# A7 Apex Locator

Input (charge the battery): DC 5V 1A  
AC Adapter: 100-240VAC 50Hz/60Hz  
Battery: NCM 14500 Rechargeable Li-ion Battery DC 3.6V/800 mAh  
Rated Power: ≤0.5W  
Measurement Voltage: AC 200mV  
Measurement Current AC: AC 100μA  
Control Unit Dimensions: L130 mm × W112 mm × H23.5 mm  
Weight: 370 gms ( including battery )

