

Superfect Zjr[®]

- Reliability
- Esthetic effect
- Biocompatibility
- Perfect fit



• HT / ST / SHT / SHTC / SHTM •



HTW[®] Superfect Zjr[®] ZIRCONIA PROFESSIONAL EXPERT

Physical parameters:


Color	white
Aesthetic	high translucency
Sintered density	≥6.0g/cm ³
Bending strength	≥1300Mpa
Fracture toughness	5Mpa ^{m^{0.5}}
Hardness	1250

HTW (high translucency white) full ceramic restoration zirconia material

[Advantage]

- Super high bending strength and excellent hardness.
- This material has comprehensive indications, especially for the long span bridge.
- The excellent masking effect can be used for covering the metal color.

[System]



[Recommended Indications]

