

# SPECTRO 70X X-RAY PERIAPICAL

**DABI ATLANTÉ**



## MOBILE COLUMN

- The Spectro 70X X-Ray Periapical Mobile Column has a base on 04 casters with super stable columns.

## WALL MOUNTED WITH SCISSOR ARMS

- The Spectro 70X X-Ray Periapical Wall mounted with scissor arms allows greater reach and use in various positions.



- With lighter and easier to move arms and a head with more precise positioning, even in the most difficult movements.
- Prepared to ensure patient and operator safety, with mechanisms that reduce the amount of exposure to the patient.
- Tube (ampoule), with focal point of 0.7 x 0.7mm, filtration with aluminum greater than 2.61 mm and winding totally immersed in special oil.
- Rated voltage 70kVp and tube current 7.0 mA.



- 5-meter manual trigger, duly tested complying with the requirements of RDC No. 330, of December 20, 2019 and IN No. 95, of May 27, 2021 of the Brazilian National Health Surveillance Agency.

Illustrative Photos. Dabi Atlante reserves the right to make changes in the design and technical specifications of the products with the authorization of ANVISA in order to pursue its improvement. Dental X-Ray Device: Anvisa 10069210087, 51000000249. Version Jun/2021.

**DABI ATLANTÉ**