



3D Pro Zir

- Strength: $\geq 1050\text{MPa}$
- Translucency: 57%

14 years
14 years of dental zirconia experience

25 years
25 years of zirconia production experience

30 years
30 years of Tosoh experience


100% TOSOH zirconia powder

[Recommended Indications]



[Colors]



[System]



3D Pro Zir Shade Matcher

Made of zirconia, the color choice of the restoration has never been so accurate.

Product features >>>



- **Bionic**
3D pro Zir Zirconia is developed by Aidite Technology and Peking University School of Stomatology. The material comes closer to natural teeth. Not only from the esthetic point of view, also its abrasion characteristics makes it to a dental material with biomimetic properties.
- **High Strength**
The cervical strength of 3D pro zir is up to 1050MPa which is greatly improved. Long span bridges are stable and fully meet the strength requirements of laboratories for a broad range of indications.
- **Layerless Natural Transition**
Matching the color transition of natural teeth from cervical to incisal, 3D pro zir has a smooth color gradient without visible layer.
- **Fast Delivery**
3D pro zir is the leading digital solution in the field of Zirconia materials. The excellent characteristics of the material itself, combined with the 3D pro zir Rapid Aesthetic Solution will help laboratories to finalize and deliver Zirconia restorations within 24 hours.